



ELK Biotechnology

CMTM6 Rabbit pAb
Catalog Number: ES6683
For research use only.

Overview

Product Name	CMTM6 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the C-terminal region of human CMTM6. at AA range; 100-180
Specificity	CMTM6 Polyclonal Antibody detects endogenous levels of CMTM6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CMTM6
Gene Name	CMTM6
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	25kDa
Human Gene ID	54918
Human Swiss-Prot Number	Q9NX76
Alternative Names	CMTM6; CKLFSF6; CKLF-like MARVEL transmembrane domain-containing protein 6; Chemokine-like factor superfamily member 6