



ELK Biotechnology

PKC δ Rabbit pAb
Catalog Number: ES6776
For research use only.

Overview

Product Name	PKC δ polyclonal antibody
Host species	Rabbit
Applications	IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunofluorescence:1/200-1/1000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human PKC δ around the non-phosphorylation site of Y52.
Specificity	PKC δ Polyclonal Antibody detects endogenous levels of PKC δ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PKC δ
Gene Name	PRKCD
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	5580
Human Swiss-Prot Number	Q05655
Alternative Names	PRKCD; Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta