

PKC γ Rabbit pAb Catalog Number: ES6783 For research use only.

Overview

Product Name	PKC γ polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/40000
Immunogen	Synthesized peptide derived from human PKC γ around the non- phosphorylation site of T655.
Specificity	PKC γ Polyclonal Antibody detects endogenous levels of PKC γ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ΡΚϹ γ
Gene Name	PRKCG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	78kDa
Human Gene ID	5582
Human Swiss-Prot Number	P05129
Alternative Names	PRKCG; PKCG; Protein kinase C gamma type; PKC-gamma