



ELK Biotechnology

PKC θ Rabbit pAb
Catalog Number: ES6796
For research use only.

Overview

Product Name	PKC θ polyclonal antibody
Host species	Rabbit
Applications	IHC-p,WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/5000
Immunogen	Synthesized peptide derived from human PKC θ around the non-phosphorylation site of S695.
Specificity	PKC θ Polyclonal Antibody detects endogenous levels of PKC θ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PKC θ
Gene Name	PRKCQ
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	76kDa
Human Gene ID	5588
Human Swiss-Prot Number	Q04759
Alternative Names	PRKCQ; PRKCT; Protein kinase C theta type; nPKC-theta