



ELK Biotechnology

DNA-PKCS Rabbit pAb
Catalog Number: ES6809
For research use only.

Overview

Product Name	DNA-PKCS polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from human DNA-PKCS around the non-phosphorylation site of S2056.
Specificity	DNA-PKCS Polyclonal Antibody detects endogenous levels of DNA-PKCS protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	DNA-PKCS
Gene Name	PRKDC
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	5591
Human Swiss-Prot Number	P78527
Alternative Names	PRKDC; HYRC; HYRC1; DNA-dependent protein kinase catalytic subunit; DNA-PK catalytic subunit; DNA-PKcs; DNPK1; p460