

MEK-1/2 (phospho Ser222/226) rabbit pAb

Cat No.:ES6825

For research use only

Overview

Product Name	MEK-1/2 (phospho Ser222/226) rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/20000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
5	peptide derived from human MEK1/2 around the
	phosphorylation site of Ser221. AA range:193-242
Specificity	Phospho-MEK-1/2 (S222/226) Polyclonal Antibody
. ,	detects endogenous levels of MEK-1/2 protein only
	when phosphorylated at S222/226.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Dual specificity mitogen-activated protein kinase
	kinase 1/2
Gene Name	MAP2K1/MAP2K2
Cellular localization	Cytoplasm, cytoskeleton, microtubule organizing
	center, centrosome . Cytoplasm, cytoskeleton,
	microtubule organizing center, spindle pole body .
	Cytoplasm . Nucleus . Membrane ; Peripheral
	membrane protein . Localizes at centrosomes during
	prometaphase, m
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	44kD
Human Gene ID	5604/5605
Human Swiss-Prot Number	Q02750/P36507



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Alternative Names

Background

MAP2K1; MEK1; PRKMK1; Dual specificity mitogen-activated protein kinase kinase 1; MAP kinase kinase 1; MAPKK 1; MKK1; ERK activator kinase 1; MAPK/ERK kinase 1; MEK 1; MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity mitogen-activated protein k The protein encoded by this gene is a member of the dual specificity protein kinase family, which acts as a mitogen-activated protein (MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein kinase lies upstream of MAP kinases and stimulates the enzymatic activity of MAP kinases upon wide variety of extra- and intracellular signals. As an essential component of MAP kinase signal transduction pathway, this kinase is involved in many cellular processes such as proliferation, differentiation, transcription regulation and development. [provided by RefSeq, Jul 2008],

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using MEK1/2 (Phospho-Ser221) Antibody. The picture on the right is blocked with the phospho peptide.





- Western blot analysis of lysates from K562 cells treated
- -- 1 with serum 20% 15' and Jurkat cells treated with EGF,
- -- 8: using MEK1/2 (Phospho-Ser221) Antibody. The lane on the right is blocked with the phospho peptide.





23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C