



ELK Biotechnology

Bax Rabbit pAb
Catalog Number: ES6948
For research use only.

Overview

Product Name	Bax polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from the N-terminal region of human Bax. at AA range; 10-90
Specificity	Bax Polyclonal Antibody detects endogenous levels of Bax protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Bax
Gene Name	BAX
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	581
Human Swiss-Prot Number	Q07815
Alternative Names	