



ELK Biotechnology

Paxillin Rabbit pAb
Catalog Number: ES6956
For research use only.

Overview

Product Name	Paxillin polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human Paxillin around the non-phosphorylation site of Y88.
Specificity	Paxillin Polyclonal Antibody detects endogenous levels of Paxillin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Paxillin
Gene Name	PXN
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	68kDa
Human Gene ID	5829
Human Swiss-Prot Number	P49023
Alternative Names	PXN; Paxillin