



ELK Biotechnology

Rad17 Rabbit pAb
Catalog Number: ES6967
For research use only.

Overview

Product Name	Rad17 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human Rad17 around the non-phosphorylation site of S646.
Specificity	Rad17 Polyclonal Antibody detects endogenous levels of Rad17 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Rad17
Gene Name	RAD17
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	77kDa
Human Gene ID	5884
Human Swiss-Prot Number	O75943
Alternative Names	RAD17; R24L; Cell cycle checkpoint protein RAD17; hRad17; RF-C/activator 1 homolog