



ELK Biotechnology

Rb (phospho Ser612) Rabbit pAb

Catalog Number: ES6997

For research use only.

Overview

Product Name	Rb (phospho Ser612) polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Immunohistochemistry:1/100-1/300 ELISA:1/5000
Immunogen	Synthesized peptide derived from human Rb around the phosphorylation site of S612.
Specificity	Phospho-Rb (S612) Polyclonal Antibody detects endogenous levels of Rb protein only when phosphorylated at S612.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Rb (phospho Ser612)
Gene Name	RB1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	5925
Human Swiss-Prot Number	P06400
Alternative Names	RB1; Retinoblastoma-associated protein; p105-Rb; pRb; Rb; pp110