

Rb Rabbit pAb Catalog Number: ES7001 For research use only.

Overview

Product Name Rb polyclonal antibody

Host species Rabbit

Applications IHC-p,IF,ELISA

Species Cross-Reactivity Human, Mouse, Rat

Recommended dilutions Immunohistochemistry:1/100-1/300

Immunofluorescence:1/200-1/1000

ELISA:1/5000

Immunogen Synthesized peptide derived from human Rb around the non-

phosphorylation site of T821.

Specificity Rb Polyclonal Antibody detects endogenous levels of Rb protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid

repeated freeze / thaw cycles.

Protein Name Rb

Gene Name RB1

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Clonality Polyclonal

Concentration 1 mg/ml

Observed band kDa

Human Gene ID 5925

Human Swiss-Prot Number P06400

Alternative Names RB1; Retinoblastoma-associated protein; p105-Rb; pRb; Rb; pp110