

Rsk-1 (phospho Thr359/S363) rabbit pAb

Cat No.: ES7081

For research use only

Overview

Product Name Rsk-1 (phospho Thr359/S363) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA **Species Cross-Reactivity** Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human p90 RSK around the phosphorylation site of Thr359 and Ser363. AA

range:331-380

Specificity Phospho-Rsk-1 (T359/S363) Polyclonal Antibody

detects endogenous levels of Rsk-1 protein only

when phosphorylated at T359/S363.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ribosomal protein S6 kinase alpha-1

Gene Name RPS6KA1

Cellular localization Nucleus. Cytoplasm.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 83kD
Human Gene ID 6195
Human Swiss-Prot Number Q15418

numan swiss-Prot Number Q15416

Alternative Names RPS6KA1; MAPKAPK1A; RSK1; Ribosomal protein S6

kinase alpha-1; S6K-alpha-1; 90 kDa ribosomal protein S6 kinase 1; p90-RSK 1; p90RSK1; p90S6K;

MAP kinase-activated protein kinase 1a;

MAPK-activated protein kinase 1a; MAPKAP kinase




+86-27-59760950 ELKbio@ELKbiotech.com




Background

1a; MAPKAP

ribosomal protein S6 kinase A1(RPS6KA1) Homo sapiens This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 nonidentical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signalling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],




Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using p90 RSK (Phospho-Thr359+Ser363) Antibody



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using p90 RSK (Phospho-Thr359+Ser363) Antibody. The picture on the right is blocked with the phospho peptide.







Western blot analysis of lysates from 293 cells treated with PMA 125ng/ml 30', using p90 RSK (Phospho-Thr359+Ser363) Antibody. The lane on the right is blocked with the phospho peptide.



+86-27-59760950