

Rsk-1 (phospho Thr573) rabbit pAb

Cat No.: ES7082

For research use only

Overview

Specificity

Product Name Rsk-1 (phospho Thr573) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications. The antiserum was produced against synthesized

Immunogen The antiserum was produced against synthesized

peptide derived from human p90 RSK around the phosphorylation site of Thr573. AA range:539-588 Phospho-Rsk-1 (T573) Polyclonal Antibody detects

endogenous levels of Rsk-1 protein only when

phosphorylated at T573.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ribosomal protein S6 kinase alpha-1

Gene Name RPS6KA1

Cellular localization Nucleus. Cytoplasm.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 95kD
Human Gene ID 6195
Human Swiss-Prot Number Q15418

Alternative Names RPS6KA1; MAPKAPK1A; RSK1; Ribosomal protein S6

kinase alpha-1; S6K-alpha-1; 90 kDa ribosomal protein S6 kinase 1; p90-RSK 1; p90RSK1; p90S6K;

MAP kinase-activated protein kinase 1a;

MAPK-activated protein kinase 1a; MAPKAP kinase



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




Background

1a; MAPKAP

ribosomal protein S6 kinase A1(RPS6KA1) Homo sapiens This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 nonidentical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signalling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],

Immunohistochemical analysis of paraffin-embedded Human colon cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absorb



0.450
0.400
0.350
0.300
0.250
0.200
0.150
0.100
0.050
0.000

phosphopeptide non-phosphopeptide


© OD 450nm Reading

Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using p90 RSK (Phospho-Thr573) Antibody






Immunofluorescence analysis of COS7 cells, using p90 RSK (Phospho-Thr573) Antibody. The picture on the right is blocked with the phospho peptide.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using p90 RSK (Phospho-Thr573) Antibody. The picture on the right is blocked with the phospho peptide.



+86-27-59760950

