



ELK Biotechnology

Sp1 Rabbit pAb
Catalog Number: ES7257
For research use only.

Overview

Product Name	Sp1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/20000
Immunogen	Synthesized peptide derived from human Sp1 around the non-phosphorylation site of T453.
Specificity	Sp1 Polyclonal Antibody detects endogenous levels of Sp1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Sp1
Gene Name	SP1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	90kDa
Human Gene ID	6667
Human Swiss-Prot Number	P08047
Alternative Names	SP1; TSFP1; Transcription factor Sp1