



ELK Biotechnology

Stat3 Rabbit pAb
Catalog Number: ES7296
For research use only.

Overview

Product Name	Stat3 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat,Monkey
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/10000
Immunogen	Synthesized peptide derived from the C-terminal region of human Stat3. at AA range; 630-710
Specificity	Stat3 Polyclonal Antibody detects endogenous levels of Stat3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Stat3
Gene Name	STAT3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	85kDa
Human Gene ID	6774
Human Swiss-Prot Number	P40763
Alternative Names	STAT3; APRF; Signal transducer and activator of transcription 3; Acute-phase response factor