

LKB1 rabbit pAb

Cat No.:ES7312

For research use only

Overview

Product Name	LKB1 rabbit pAb	
Host species	Rabbit	
Applications	WB;IHC;IF;ELISA	
Species Cross-Reactivity	Human;Mouse	
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300.	
	Immunofluorescence: 1/200 - 1/1000. ELISA:	
	1/20000. Not yet tested in other applications.	
Immunogen	The antiserum was produced against synthesized	
	peptide derived from human LKB1. AA	
	range:384-433	
Specificity	LKB1 Polyclonal Antibody detects endogenous levels	
	of LKB1 protein.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	Serine/threonine-protein kinase STK11	
Gene Name	STK11	
Cellular localization	Nucleus. Cytoplasm. Membrane . Mitochondrion. A	
	small fraction localizes at membranes (By similarity).	
	Relocates to the cytoplasm when bound to STRAD	
	(STRADA or STRADB) and CAB39/MO25	
	(CAB39/MO25alpha or CAB39L/MO25beta).	
	Translocates to the mitochondrion	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band		
Human Gene ID	6794	
Human Swiss-Prot Number	Q15831	
Alternative Names	STK11; LKB1; PJS; Serine/threonine-protein kinase	
	STK11; Liver kinase B1; LKB1; hLKB1; Renal	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

colo 293T

138-

70-

40---35---25---


15---

carcinoma antigen NY-REN-19

This gene, which encodes a member of the serine/threonine kinase family, regulates cell polarity and functions as a tumor suppressor. Mutations in this gene have been associated with Peutz-Jeghers syndrome, an autosomal dominant disorder characterized by the growth of polyps in the gastrointestinal tract, pigmented macules on the skin and mouth, and other neoplasms. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized. [provided by RefSeq, Jul 2008],

Western blot analysis of colo 293T lysis using LKB1 antibody. Antibody was diluted at 1:500

Immunofluorescence analysis of NIH/3T3 cells, using LKB1 Antibody. The picture on the right is blocked with the synthesized peptide.





+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using LKB1 Antibody. The picture on the right is blocked with the synthesized peptide.





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C