



ELK Biotechnology

BST-2 Rabbit pAb
Catalog Number: ES7317
For research use only.

Overview

Product Name	BST-2 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human BST-2. at AA range; 50-130
Specificity	BST-2 Polyclonal Antibody detects endogenous levels of BST-2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	BST-2
Gene Name	BST2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	22kDa
Human Gene ID	684
Human Swiss-Prot Number	Q10589
Alternative Names	BST2; Bone marrow stromal antigen 2; BST-2; HM1.24 antigen; Tetherin; CD antigen CD317