



# ELK Biotechnology

**p53 Rabbit pAb**  
**Catalog Number: ES7433**  
**For research use only.**

## Overview

---

<b>Product Name</b>	p53 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	WB,IHC-p,IF,ELISA
<b>Species Cross-Reactivity</b>	Human,Mouse,Rat,Monkey
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>Immunofluorescence:1/200-1/1000</b> <b>ELISA:1/20000</b>
<b>Immunogen</b>	Synthesized peptide derived from human p53 around the non-acetylation site of K386. at AA range: 320-400
<b>Specificity</b>	p53 Polyclonal Antibody detects endogenous levels of p53 protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	p53
<b>Gene Name</b>	TP53
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>53kDa</b>
<b>Human Gene ID</b>	7157
<b>Human Swiss-Prot Number</b>	P04637
<b>Alternative Names</b>	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53