



ELK Biotechnology

p73 Rabbit pAb
Catalog Number: ES7440
For research use only.

Overview

Product Name	p73 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from the Internal region of human p73. at AA range; 240-320
Specificity	p73 Polyclonal Antibody detects endogenous levels of p73 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	p73
Gene Name	TP73
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	70kDa
Human Gene ID	7161
Human Swiss-Prot Number	O15350
Alternative Names	TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein