



# ELK Biotechnology

**Ox40L Rabbit pAb**  
**Catalog Number: ES7466**  
**For research use only.**

## Overview

---

<b>Product Name</b>	Ox40L polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>IF,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human</b>
<b>Recommended dilutions</b>	<b>Immunofluorescence:1/200-1/1000</b> <b>ELISA:1/20000</b>
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of human Ox40L. at AA range; 50-130
<b>Specificity</b>	Ox40L Polyclonal Antibody detects endogenous levels of Ox40L protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	Ox40L
<b>Gene Name</b>	TNFSF4
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>kDa</b>
<b>Human Gene ID</b>	7292
<b>Human Swiss-Prot Number</b>	P23510
<b>Alternative Names</b>	TNFSF4; TXGP1; Tumor necrosis factor ligand superfamily member 4; Glycoprotein Gp34; OX40 ligand; OX40L; TAX transcriptionally-activated glycoprotein 1; CD antigen CD252