



ELK Biotechnology

ZAP-70 Rabbit pAb
Catalog Number: ES7535
For research use only.

Overview

Product Name	ZAP-70 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from human ZAP-70 around the non-phosphorylation site of Y292.
Specificity	ZAP-70 Polyclonal Antibody detects endogenous levels of ZAP-70 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ZAP-70
Gene Name	ZAP70
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	70kDa
Human Gene ID	7535
Human Swiss-Prot Number	P43403
Alternative Names	ZAP70; SRK; Tyrosine-protein kinase ZAP-70; 70 kDa zeta-chain associated protein; Syk-related tyrosine kinase