



ELK Biotechnology

Frizzled-3 Rabbit pAb
Catalog Number: ES7593
For research use only.

Overview

| | |
|--------------------------|---|
| Product Name | Frizzled-3 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IF,ELISA |
| Species Cross-Reactivity | Human,Mouse |
| Recommended dilutions | Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/40000 |
| Immunogen | Synthesized peptide derived from the Internal region of human Frizzled-3. at AA range; 460-540 |
| Specificity | Frizzled-3 Polyclonal Antibody detects endogenous levels of Frizzled-3 protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | Frizzled-3 |
| Gene Name | FZD3 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 76kDa |
| Human Gene ID | 7976 |
| Human Swiss-Prot Number | Q9NPG1 |
| Alternative Names | FZD3; Frizzled-3; Fz-3; hFz3 |