

Olfactory receptor 2AE1 rabbit pAb

Cat No.: ES7634

For research use only

Overview

Product Name Olfactory receptor 2AE1 rabbit pAb

Host species Rabbit
Applications IF;ELISA

Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions Immunofluorescence: 1/200 - 1/1000. ELISA:

1/10000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human OR2AE1. AA

range:231-280

Specificity Olfactory receptor 2AE1 Polyclonal Antibody detects

endogenous levels of Olfactory receptor 2AE1

protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Olfactory receptor 2AE1

Gene Name OR2AE1

Cellular localizationCell membrane; Multi-pass membrane protein.PurificationThe antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 81392 Human Swiss-Prot Number Q8NHA4

Alternative Names OR2AE1; OR2AE2; Olfactory receptor 2AE1;

Olfactory receptor 2AE2

Background Olfactory receptors interact with odorant molecules

in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from





single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],

Immunofluorescence analysis of A549 cells, using OR2AE1 Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950

