



# ELK Biotechnology

## Olfactory receptor 2G3 Rabbit pAb

Catalog Number: ES7637

For research use only.

### Overview

---

Product Name	Olfactory receptor 2G3 polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>Immunofluorescence:1/200-1/1000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from the Internal region of human Olfactory receptor 2G3. at AA range; 30-110
Specificity	Olfactory receptor 2G3 Polyclonal Antibody detects endogenous levels of Olfactory receptor 2G3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Olfactory receptor 2G3
Gene Name	OR2G3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>35kDa</b>
Human Gene ID	81469
Human Swiss-Prot Number	Q8NGZ4
Alternative Names	OR2G3; Olfactory receptor 2G3; Olfactory receptor OR1-33