



ELK Biotechnology

ELL Rabbit pAb
Catalog Number: ES7650
For research use only.

Overview

Product Name	ELL polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from the Internal region of human ELL. at AA range; 60-140
Specificity	ELL Polyclonal Antibody detects endogenous levels of ELL protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ELL
Gene Name	ELL
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	65kDa
Human Gene ID	8178
Human Swiss-Prot Number	P55199
Alternative Names	ELL; C19orf17; RNA polymerase II elongation factor ELL; Eleven-nineteen lysine-rich leukemia protein