



# ELK Biotechnology

**NCoA-3 Rabbit pAb**  
**Catalog Number: ES7656**  
**For research use only.**

## Overview

---

<b>Product Name</b>	NCoA-3 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>IHC-p,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse</b>
<b>Recommended dilutions</b>	<b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/40000</b>
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of human NCoA-3.
<b>Specificity</b>	NCoA-3 Polyclonal Antibody detects endogenous levels of NCoA-3 protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	NCoA-3
<b>Gene Name</b>	NCOA3
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>kDa</b>
<b>Human Gene ID</b>	8202
<b>Human Swiss-Prot Number</b>	Q9Y6Q9
<b>Alternative Names</b>	NCOA3; AIB1; BHLHE42; RAC3; TRAM1; Nuclear receptor coactivator 3; NCoA-3; ACTR; Amplified in breast cancer 1 protein; AIB-1; CBP-interacting protein; pCIP; Class E basic helix-loop-helix protein 42; bHLHe42; Receptor-associated coactivator