



ELK Biotechnology

Frizzled-1 Rabbit pAb
Catalog Number: ES7664
For research use only.

Overview

Product Name	Frizzled-1 polyclonal antibody
Host species	Rabbit
Applications	WB,IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunofluorescence:1/200-1/1000 ELISA:1/20000
Immunogen	Synthesized peptide derived from the N-terminal region of human Frizzled-1. at AA range; 30-110
Specificity	Frizzled-1 Polyclonal Antibody detects endogenous levels of Frizzled-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Frizzled-1
Gene Name	FZD1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	71kDa
Human Gene ID	8321
Human Swiss-Prot Number	Q9UP38
Alternative Names	FZD1; Frizzled-1; Fz-1; hFz1; FzE1