



# ELK Biotechnology

**Frizzled-4 Rabbit pAb**  
**Catalog Number: ES7666**  
**For research use only.**

## Overview

---

Product Name	Frizzled-4 polyclonal antibody
Host species	Rabbit
Applications	IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	<b>Immunofluorescence:1/200-1/1000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from the Internal region of human Frizzled-4. at AA range; 100-180
Specificity	Frizzled-4 Polyclonal Antibody detects endogenous levels of Frizzled-4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Frizzled-4
Gene Name	FZD4
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>kDa</b>
Human Gene ID	8322
Human Swiss-Prot Number	Q9ULV1
Alternative Names	FZD4; Frizzled-4; Fz-4; hFz4; FzE4; CD antigen CD344