



# ELK Biotechnology

**Histone H3 Rabbit pAb**  
**Catalog Number: ES7688**  
**For research use only.**

## Overview

---

<b>Product Name</b>	Histone H3 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,IHC-p,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/5000</b>
<b>Immunogen</b>	Synthesized peptide derived from human Histone H3 around the non-phosphorylation site of T3.
<b>Specificity</b>	Histone H3 Polyclonal Antibody detects endogenous levels of Histone H3 protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	Histone H3
<b>Gene Name</b>	HIST1H3A
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>22kDa</b>
<b>Human Gene ID</b>	8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
<b>Human Swiss-Prot Number</b>	P68431
<b>Alternative Names</b>	HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3