

## GPR101 Rabbit pAb Catalog Number: ES7696 For research use only.

## Overview

Product Name GPR101 polyclonal antibody

Host species Rabbit

Applications WB,IF,ELISA

Species Cross-Reactivity Human

Recommended dilutions Western Blot:1/500-1/2000

Immunofluorescence: 1/200-1/1000

ELISA:1/5000

**Immunogen** Synthesized peptide derived from the Internal region of human GPR101. at

AA rangle; 220-300

**Specificity** GPR101 Polyclonal Antibody detects endogenous levels of GPR101 protein.

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid

repeated freeze / thaw cycles.

Protein Name GPR101

Gene Name GPR101

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Clonality** Polyclonal

Concentration 1 mg/ml

Observed band 56kDa

Human Gene ID 83550

Human Swiss-Prot Number Q96P66

Alternative Names GPR101; Probable G-protein coupled receptor 101

E: elkbio@elkbiotech.com