

## NFATc2IP rabbit pAb

Cat No.: ES7785

For research use only

## Overview

Product Name NFATc2IP rabbit pAb

Host species Rabbit
Applications WB;IHC

**Species Cross-Reactivity** Human;Rat;Mouse;;Fish

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300

Immunogen The antiserum was produced against synthesized

peptide derived from human NFATC2IP. AA

range:151-200

Specificity NFATc2IP Polyclonal Antibody detects endogenous

levels of NFATc2IP protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Protein Name** NFATC2-interacting protein

Gene Name NFATC2IP

Cellular localization Nucleus . Cytoplasm . TRAF1 is associated with a

fraction of NFATC2IP in the cytoplasm and prevents

its translocation to the nucleus. .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 48kD
Human Gene ID 84901
Human Swiss-Prot Number Q8NCF5

Alternative Names NFATC2IP; NIP45; NFATC2-interacting protein; 45

kDa NF-AT-interacting protein; 45 kDa

NFAT-interacting protein; Nuclear factor of activated

T-cells; cytoplasmic 2-interacting protein

**Background** function:Plays a role in the inducible expression of

cytokine genes in T-cells, especially by increasing NFAT-driven IL-4 production., PTM: Methylation at



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbio



the N-terminus by PRMT1 modulates interaction with the NFAT complex and results in augmented cytokine production.,similarity:Contains 1 ubiquitin-like domain.,subcellular location:TRAF1 is associated with a fraction of NFATC2IP in the cytoplasm and prevents its translocation to the nucleus.,subunit:Interacts with NFATC2, TRAF1, TRAF2 and PRMT1.,



+86-27-59760950

Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

