

CHST1 rabbit pAb

Cat No.: ES7811

For research use only

Overview

Product Name CHST1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not

yet tested in other applications.

Immunogen Synthesized peptide derived from the Internal region

of human CHST1.

Specificity CHST1 Polyclonal Antibody detects endogenous

levels of CHST1 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Carbohydrate sulfotransferase 1

Gene Name CHST1

Cellular localization Golgi apparatus membrane; Single-pass type II

membrane protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 49kD
Human Gene ID 8534
Human Swiss-Prot Number 043916

Alternative Names CHST1; Carbohydrate sulfotransferase 1;

Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 1; GST-1; Keratan sulfate Gal-6

 $sulfot ransferase; \, KS6ST; \, KSGal6ST; \, KSST$

Background This locus encodes a member of the keratin

sulfotransferase family of proteins. The encoded enzyme catalyzes the sulfation of the proteoglycan

keratin. [provided by RefSeq, Aug 2011],



+86-27-59760950 ELKbio@ELKbiotech.com