

## RUNX1 Rabbit pAb Catalog Number: ES7836 For research use only.

## Overview

Product Name RUNX1 polyclonal antibody

Host species Rabbit

Applications WB,ELISA

Species Cross-Reactivity Human, Mouse, Rat

Recommended dilutions Western Blot:1/500-1/2000

ELISA:1/5000

Immunogen Synthesized peptide derived from human RUNX1 around the non-

phosphorylation site of S435.

**Specificity** RUNX1 Polyclonal Antibody detects endogenous levels of RUNX1 protein.

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid

repeated freeze / thaw cycles.

Protein Name RUNX1

Gene Name RUNX1

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Clonality** Polyclonal

Concentration 1 mg/ml

Observed band 53kDa

Human Gene ID 861

Human Swiss-Prot Number Q01196

Alternative Names RUNX1; AML1; CBFA2; Runt-related transcription factor 1; Acute myeloid

leukemia 1 protein; Core-binding factor subunit alpha-2; CBF-alpha-2; Oncogene AML-1; Polyomavirus enhancer-binding protein 2 alpha B

subunit; PEA2-alpha B; PEBP2-alpha

T:86-27-59760950