



ELK Biotechnology

PEA-15 Rabbit pAb
Catalog Number: ES7858
For research use only.

Overview

Product Name	PEA-15 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/5000
Immunogen	Synthesized peptide derived from human PEA-15 around the non-phosphorylation site of S104.
Specificity	PEA-15 Polyclonal Antibody detects endogenous levels of PEA-15 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PEA-15
Gene Name	PEA15
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	15kDa
Human Gene ID	8682
Human Swiss-Prot Number	Q15121
Alternative Names	PEA15; Astrocytic phosphoprotein PEA-15; 15 kDa phosphoprotein enriched in astrocytes; Phosphoprotein enriched in diabetes; PED