



ELK Biotechnology

IL-18R Rabbit pAb
Catalog Number: ES7880
For research use only.

Overview

| | |
|---------------------------------|--|
| Product Name | IL-18R polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,ELISA |
| Species Cross-Reactivity | Human |
| Recommended dilutions | Western Blot:1/500-1/2000 ELISA:1/5000 |
| Immunogen | Synthesized peptide derived from the Internal region of human IL-18R. at AA range; 200-280 |
| Specificity | IL-18R Polyclonal Antibody detects endogenous levels of IL-18R protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | IL-18R |
| Gene Name | IL18R1 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 65kDa |
| Human Gene ID | 8809 |
| Human Swiss-Prot Number | Q13478 |
| Alternative Names | IL18R1; IL1RRP; Interleukin-18 receptor 1; IL-18R-1; IL-18R1; CD218 antigen-like family member A; CDw218a; IL1 receptor-related protein; IL-1Rrp; IL1R-rp; CD antigen CD218a |