



ELK Biotechnology

Neuropilin-2 Rabbit pAb Catalog Number: ES7884 For research use only.

Overview

Product Name	Neuropilin-2 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/40000
Immunogen	Synthesized peptide derived from the Internal region of human Neuropilin-2.
Specificity	Neuropilin-2 Polyclonal Antibody detects endogenous levels of Neuropilin-2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Neuropilin-2
Gene Name	NRP2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	120kDa
Human Gene ID	8828
Human Swiss-Prot Number	O60462
Alternative Names	NRP2; VEGF165R2; Neuropilin-2; Vascular endothelial cell growth factor 165 receptor 2