



ELK Biotechnology

Claudin-6 Rabbit pAb Catalog Number: ES7951 For research use only.

Overview

Product Name	Claudin-6 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Immunohistochemistry:1/100-1/300 ELISA:1/5000
Immunogen	Synthesized peptide derived from human Claudin-6 around the non-phosphorylation site of Y219.
Specificity	Claudin-6 Polyclonal Antibody detects endogenous levels of Claudin-6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Claudin-6
Gene Name	CLDN6
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	9074
Human Swiss-Prot Number	P56747
Alternative Names	CLDN6; Claudin-6; Skullin