



ELK Biotechnology

MAPKAPK-2 Rabbit pAb
Catalog Number: ES8006
For research use only.

Overview

| | |
|---------------------------------|--|
| Product Name | MAPKAPK-2 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000 |
| Immunogen | Synthesized peptide derived from human MAPKAPK-2 around the non-phosphorylation site of T222. |
| Specificity | MAPKAPK-2 Polyclonal Antibody detects endogenous levels of MAPKAPK-2 protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | MAPKAPK-2 |
| Gene Name | MAPKAPK2 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 46kDa |
| Human Gene ID | 9261 |
| Human Swiss-Prot Number | P49137 |
| Alternative Names | MAPKAPK2; MAP kinase-activated protein kinase 2; MAPK-activated protein kinase 2; MAPKAP kinase 2; MAPKAP-K2; MAPKAPK-2; MK-2; MK2 |