

## Gab 2 (phospho Tyr643) rabbit pAb

**Cat No.:ES8108** 

For research use only

## Overview

Product Name Gab 2 (phospho Tyr643) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

**Recommended dilutions** Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications. The antiserum was produced against synthesized

Immunogen The antiserum was produced against synthesized

peptide derived from human Gab2 around the phosphorylation site of Tyr643. AA range:609-658

**Specificity** Phospho-Gab 2 (Y643) Polyclonal Antibody detects

endogenous levels of Gab 2 protein only when

phosphorylated at Y643.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Protein Name** GRB2-associated-binding protein 2

Gene Name GAB2

**Cellular localization** Cytoplasm . Cell membrane .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 74kD
Human Gene ID 9846
Human Swiss-Prot Number Q9UQC2

Alternative Names GAB2; KIAA0571; GRB2-associated-binding protein

2; GRB2-associated binder 2; Growth factor receptor

bound protein 2-associated protein 2; pp100 GRB2 associated binding protein 2(GAB2) Homo

sapiens This gene is a member of the

GRB2-associated binding protein (GAB) gene family.




**Background** 


+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com




These proteins contain pleckstrin homology (PH) domain, and bind SHP2 tyrosine phosphatase and GRB2 adapter protein. They act as adapters for transmitting various signals in response to stimuli through cytokine and growth factor receptors, and T- and B-cell antigen receptors. The protein encoded by this gene is the principal activator of phosphatidylinositol-3 kinase in response to activation of the high affinity IgE receptor. Two alternatively spliced transcripts encoding different isoforms have been described for this gene. [provided by RefSeq, Nov 2009],

Immunohistochemistry analysis of paraffin-embedded human brain, using Gab2 (Phospho-Tyr643) Antibody. The picture on the right is blocked with the phospho peptide.







Western blot analysis of lysates from Jurkat cells treated with IFN 2500U/ML 30', using Gab2 (Phospho-Tyr643) Antibody. The lane on the right is blocked with the phospho peptide.

