



ELK Biotechnology

beta 1 Sodium Potassium ATPase Rabbit pAb

Catalog Number: ES8173

For research use only.

Overview

Product Name	beta 1 Sodium Potassium ATPase polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/2000
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous beta 1 Sodium Potassium ATPase protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	beta 1 Sodium Potassium ATPase
Gene Name	ATP1B1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	45-50kDa
Human Gene ID	481
Human Swiss-Prot Number	P05026
Alternative Names	