



# ELK Biotechnology

**$\beta$ -Catenin Rabbit pAb**  
**Catalog Number: ES8174**  
**For research use only.**

## Overview

---

Product Name	$\beta$ -Catenin polyclonal antibody
Host species	Rabbit
Applications	<b>WB</b>
Species Cross-Reactivity	<b>Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/2000-1/4000</b>
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous $\beta$ -Catenin protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	$\beta$ -Catenin
Gene Name	CTNNB1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>92kDa</b>
Human Gene ID	1499
Human Swiss-Prot Number	P35222
Alternative Names	CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin