



ELK Biotechnology

PEPC Rabbit pAb
Catalog Number: ES8175
For research use only.

Overview

Product Name	PEPC polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Plants
Recommended dilutions	Western Blot:1/2000-1/5000
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous PEPC protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PEPC
Gene Name	PGC
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	105-109kDa
Human Gene ID	5225
Human Swiss-Prot Number	P20142
Alternative Names	PGC; Gastricsin; Pepsinogen C