



ELK Biotechnology

E-cadherin Rabbit pAb
Catalog Number: ES8181
For research use only.

Overview

Product Name	E-cadherin polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/1000 Immunohistochemistry:1/500-1/1000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous E-cadherin protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	E-cadherin
Gene Name	CDH1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	135kDa
Human Gene ID	999
Human Swiss-Prot Number	P12830
Alternative Names	