



# ELK Biotechnology

**14-3-3  $\epsilon$  Rabbit pAb**  
**Catalog Number: ES8183**  
**For research use only.**

## Overview

---

Product Name	14-3-3 $\epsilon$ polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Mouse,Rat
Recommended dilutions	<b>Western Blot:1/1000</b> <b>Immunohistochemistry:1/200-1/500</b>
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous 14-3-3 $\epsilon$ protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	14-3-3 $\epsilon$
Gene Name	YWHAE
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>28kDa</b>
Human Gene ID	7531
Human Swiss-Prot Number	P62258
Alternative Names	