



# ELK Biotechnology

**HSP60 Rabbit pAb**  
**Catalog Number: ES8185**  
**For research use only.**

## Overview

---

Product Name	HSP60 polyclonal antibody
Host species	Rabbit
Applications	<b>WB,IHC-p</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/1000-1/2000</b> <b>Immunohistochemistry:1/200-1/500</b>
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous HSP60 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HSP60
Gene Name	HSPD1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>60kDa</b>
Human Gene ID	3329
Human Swiss-Prot Number	P10809
Alternative Names	