



ELK Biotechnology

HSC 70 Rabbit pAb
Catalog Number: ES8186
For research use only.

Overview

Product Name	HSC 70 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human HSC 70.
Specificity	HSC 70 Polyclonal Antibody detects endogenous levels of HSC 70 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HSC 70
Gene Name	HSPA8
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	70-72kDa
Human Gene ID	3312
Human Swiss-Prot Number	P11142
Alternative Names	HSPA8; HSC70; HSP73; HSPA10; Heat shock cognate 71 kDa protein; Heat shock 70 kDa protein 8