



# ELK Biotechnology

**HIF-1 $\beta$  Rabbit pAb**  
**Catalog Number: ES8188**  
**For research use only.**

## Overview

---

<b>Product Name</b>	HIF-1 $\beta$ polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,IHC-p</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/1000-1/2000</b> <b>Immunohistochemistry:1/200-1/500</b>
<b>Immunogen</b>	Recombinant Protein
<b>Specificity</b>	The antibody detects endogenous HIF-1 $\beta$ protein.
<b>Formulation</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	HIF-1 $\beta$
<b>Gene Name</b>	ARNT
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>87kDa</b>
<b>Human Gene ID</b>	405
<b>Human Swiss-Prot Number</b>	P27540
<b>Alternative Names</b>	