



ELK Biotechnology

ERK 3 Rabbit pAb
Catalog Number: ES8192
For research use only.

Overview

Product Name	ERK 3 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/1000 Immunohistochemistry:1/200-1/500
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous ERK 3 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ERK 3
Gene Name	MAPK6
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	83-105kDa
Human Gene ID	5597
Human Swiss-Prot Number	Q16659
Alternative Names	MAPK6; ERK3; PRKM6; Mitogen-activated protein kinase 6; MAP kinase 6; MAPK 6; Extracellular signal-regulated kinase 3; ERK-3; MAP kinase isoform p97; p97-MAPK