



# ELK Biotechnology

**Histone H4 Rabbit pAb**  
**Catalog Number: ES8193**  
**For research use only.**

## Overview

---

<b>Product Name</b>	Histone H4 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b>
<b>Immunogen</b>	Synthetic Peptide
<b>Specificity</b>	The antibody detects endogenous Histone H4 protein.
<b>Formulation</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	Histone H4
<b>Gene Name</b>	HIST1H4A; HIST1H4B; HIST1H4C; HIST1H4D; HIST1H4E; HIST1H4F; HIST1H4H; HIST1H4I; HIST1H4J; HIST1H4K; HIST1H4L; HIST2H4A; HIST2H4B; HIST4H4
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>14kDa</b>
<b>Human Gene ID</b>	8294
<b>Human Swiss-Prot Number</b>	P62805
<b>Alternative Names</b>	