



# ELK Biotechnology

**GSK3 $\beta$  Rabbit pAb**  
**Catalog Number: ES8198**  
**For research use only.**

## Overview

---

|                          |   |
|--------------------------|---|
| Product Name             | GSK3 $\beta$ polyclonal antibody  |
| Host species             | Rabbit  |
| Applications             | WB,IHC-p  |
| Species Cross-Reactivity | Mouse,Rat   |
| Recommended dilutions    | <b>Western Blot:1/1000-1/2000</b><br><b>Immunohistochemistry:1/200-1/500</b>  |
| Immunogen                | Recombinant Protein   |
| Specificity              | The antibody detects endogenous GSK3 $\beta$ protein.   |
| Formulation              | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.  |
| Storage                  | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.               |
| Protein Name             | GSK3 $\beta$  |
| Gene Name                | GSK3B   |
| Purification             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality                | Polyclonal  |
| Concentration            | 1 mg/ml   |
| Observed band            | <b>46kDa</b>  |
| Human Gene ID            | 2932  |
| Human Swiss-Prot Number  | P49841  |
| Alternative Names        |   |